

# session\_demonstration\_script

November 10, 2021

## 1 Example code for using session.Session

**Note:** This notebook covers several relevant methods of the `Session` and `Stim` objects, detailing some of their arguments, as well. For more details, take a look at the docstring associated with a method of interest.

### Import notes:

- These packages should be present if installing the conda environment from `osca.yml`.
- `util` is a [Github repo](#) of mine, and the correct branch `osca_mult` is automatically installed from `osca.yml`. Errors internal to the codebase involving `util` code and occurring *after* new changes have been pulled from the `OpenScope_CA_Analysis` repo *may* be due to an update of the `osca_mult` branch of `util`. Though I will try to avoid this, consider updating the utility under those circumstances.

### 1.1 Set paths to main data directory and the mouse dataframe

If you wish to use the same formatting style (and logging format) as I do:

### 1.2 1. Basics of initializing a Session object

After creating the session, you must run `self.extract_info()`. This wasn't amalgamated into the `__init__` to reduce the amount of information needed to just create a session object.

#### 1.2.1 Loading ROI/running/pupil info

You can load this information when you call `self.extract_info()` or manually later by calling `self.load_roi_info()`, `self.load_run_data()` and `self.load_pup_data()`.

Loading stimulus and alignment info...

Loading ROI trace info...

WARNING: Session 764704289: 3 noisy ROIs (mean below 0 or median above midrange) are also included in the NaN ROI attributes (but not set to NaN): 244, 298, 305.

Loading running info...

WARNING: Session 764704289: 211 dropped running frames (~0.1%) (in pre-processing).

Loading pupil info...

```

[5]:      stimulus_type stimulus_template_name unexpected gabor_frame \
0      grayscreen      grayscreen      NaN
1      gabors      gabors      0.0      A
2      gabors      gabors      0.0      B
3      gabors      gabors      0.0      C
4      gabors      gabors      0.0      D
...      ...
8839      visflow      visflow_right      0.0
8840      visflow      visflow_right      1.0
8841      visflow      visflow_right      1.0
8842      visflow      visflow_right      1.0
8843      grayscreen      grayscreen      NaN

      gabor_kappa gabor_mean_orientation gabor_number \
0      NaN      NaN      NaN
1      16.0      135.0      30.0
2      16.0      135.0      30.0
3      16.0      135.0      30.0
4      16.0      135.0      30.0
...      ...
8839      NaN      NaN      NaN
8840      NaN      NaN      NaN
8841      NaN      NaN      NaN
8842      NaN      NaN      NaN
8843      NaN      NaN      NaN

      gabor_locations_x \
0      []
1      [-0.998732530996428, -0.7988942745979938, -0.0...
2      [-0.8273358833992613, -0.32202169430120714, -0...
3      [-0.1439318404380644, -0.9639223437829889, -0...
4      [-0.1981534893873622, -0.7603480104179756, -0...
...      ...
8839      []
8840      []
8841      []
8842      []
8843      []

      gabor_locations_y \
0      []
1      [-0.936204215614872, -0.48115197167416995, -0...
2      [-0.6021449948480063, -0.6653905125829843, -0...
3      [-0.725093701321675, -0.021010443830197678, -0...
4      [-0.23152862741244445, -0.41875478323604776, -...
...      ...
8839      []

```

8840	[]
8841	[]
8842	[]
8843	[]

	gabor_sizes	...	\
0		[]	...
1	[293, 392, 392, 323, 280, 396, 316, 363, 226, ...	...	
2	[313, 319, 262, 228, 400, 210, 264, 218, 308, ...	...	
3	[396, 212, 277, 210, 390, 329, 406, 317, 358, ...	...	
4	[326, 244, 208, 212, 251, 242, 341, 299, 406, ...	...	
...		...	...
8839		[]	...
8840		[]	...
8841		[]	...
8842		[]	...
8843		[]	...

	square_proportion_flipped	start_frame_stim	stop_frame_stim	\
0	NaN	0	1800	
1	NaN	1800	1818	
2	NaN	1818	1836	
3	NaN	1836	1854	
4	NaN	1854	1872	
...	...	...	...	
8839	0.00	249960	250020	
8840	0.25	250020	250080	
8841	0.25	250080	250140	
8842	0.25	250140	250200	
8843	NaN	250200	251999	

	num_frames_stim	start_frame_twop	stop_frame_twop	num_frames_twop	\
0	1800	143	1046	903	
1	18	1046	1055	9	
2	18	1055	1064	9	
3	18	1064	1073	9	
4	18	1073	1082	9	
...	...	...	...	...	
8839	60	125552	125582	30	
8840	60	125582	125612	30	
8841	60	125612	125642	30	
8842	60	125642	125672	30	
8843	1799	125672	126575	903	

	start_time_sec	stop_time_sec	duration_sec
0	14.277090	44.301717	30.024627
1	44.301717	44.602241	0.300524

2	44.602241	44.902563	0.300322
3	44.902563	45.202768	0.300204
4	45.202768	45.503007	0.300240
...	...	...	...
8839	4183.741890	4184.742721	1.000831
8840	4184.742721	4185.743529	1.000808
8841	4185.743529	4186.744364	1.000835
8842	4186.744364	4187.745223	1.000860
8843	4187.745223	4217.728557	29.983333

[8844 rows x 24 columns]

### 1.2.2 Some information contained in the session object

Note: Stim objects (subclasses: Gabors, Visflow, Grayscr) are a separate class from Session objects. However, each can be accessed from the other using: - from Session objects: `self.stims`, `self.gabors`, `self.visflow` - from Stim objects: `self.sess`

```
number of rois: 628
mouse number: 6
mouse ID: 413663
gabor object: Gabors (stimulus from session 764704289)
2p frames per sec: 30.08
stimulus frames per sec: 59.95
```

### 1.3 2. Identifying segments of interest

From a Session's Stim, you can get a list of segments that fit a specific criterion, e.g. **U segments** (unexpected, 3rd segment).

Then, you can access the frame numbers.

**Note:** Specifying `ch_fl` (check flanks) ensures that only frame numbers whose flanks are within the recording are returned. In other words, any frame number too close to the start or end of the recording (based on `pre/post` values), will be dropped.

You can now get the **ROI/running/pupil data** corresponding to these reference frames and specified `pre/post` periods (in sec).

You can also directly obtain statistics on the data of interest

```
[11]: datatype          roi_traces
nan_rois_removed      yes
scaled                no
baseline              no
integrated            yes
smoothing             no
fluorescence          dff
general ROIs sequences
stats    None stat_mean 0.026752
```

```
error_SEM    0.000911
```

Data and statistics are returned in a hierarchical dataframe with **columns** and **indices**.

This has the advantage of allowing metadata to be stored in dummy columns, however extracting data from these dataframes can be tricky, syntactically.

```
[12]: datatype          roi_traces
      nan_rois_removed      yes
      scaled                yes
      baseline              no
      integrated            no
      smoothing             no
      fluorescence         dff
      ROIs sequences time_values
0      0          -1.000000   -0.009556
      -0.966102   -0.644810
      -0.932203   -0.214521
      -0.898305   -0.116127
      -0.864407   -0.318214
...
643  95           0.864407    0.050568
      0.898305    0.445153
      0.932203    0.108850
      0.966102    0.116475
      1.000000    0.213779
```

```
[3617280 rows x 1 columns]
```


To **extract a numpy array** with the correct dimensions from a hierarchical dataframe, you can use the following utility.

Here, each index level, then column level becomes an axis, **i.e. ROIs x sequences x time\_values** (In this case, `squeeze_cols` is set to `True` to prevent each dummy column from becoming an axis.)

```
ROI data shape: 628 ROIs x 96 sequences x 60 time values
```


You can also retrieve the time stamps for each frame.

Finally, we can plot each ROIs mean activity across sequences, as well as a mean across ROIs.



### 1.3.1 The same steps apply for Visflow

```
[16]: <module 'sess_util.sess_data_util' from '../sess_util/sess_data_util.py'>
```



## 1.4 3. Additional tips on indexing a hierarchical dataframe

```
[18]: scaled          yes
      baseline        no
      integrated      no
      smoothing       no
      fluorescence    dff
      ROIs sequences time_values
0    0          -1.000000    0.045433
      -0.966102    -0.303140
      -0.932203    -0.301464
```

		-0.898305	-0.219794
		-0.864407	0.065587
...			...
643	4	0.864407	-0.513739
		0.898305	-0.143983
		0.932203	0.000693
		0.966102	0.140502
		1.000000	0.283351

[113040 rows x 1 columns]

#### 1.4.1 4. Retrieving several Session objects, based on criteria

This function keeps track of which Sessions or Mice must be left out (e.g., due to a problem with the session data or the mouse didn't see the stimulus of interest <- the latter only comes up with pilot data).

You can now retrieve the mouse number, session number and ID that fit specific the criteria,

e.g., **session number 1, 2 or 3, production, dendritic plane**

WARNING: Sorted and unique will be set to False as multiple labels are requested.

```

mouse 6: 764704289 (session 1)
mouse 6: 765193831 (session 2)
mouse 6: 766502238 (session 3)
mouse 8: 777914830 (session 1)
mouse 8: 778864809 (session 2)
mouse 8: 779650018 (session 3)
mouse 9: 826187862 (session 1)
mouse 9: 826773996 (session 2)
mouse 9: 827833392 (session 3)
mouse 10: 826338612 (session 1)
mouse 10: 826819032 (session 2)
mouse 10: 828816509 (session 3)
mouse 11: 823453391 (session 1)
mouse 11: 824434038 (session 2)
mouse 11: 825180479 (session 3)

```

You can now **initialize the Sessions** using this function which does the additional extraction steps automatically.

Creating session 764704289...

Loading stimulus and alignment info...

Loading ROI trace info...

WARNING: Session 764704289: 3 noisy ROIs (mean below 0 or median above midrange) are also included in the NaN ROI attributes (but not set to NaN): 244, 298, 305.

Loading running info...

WARNING: Session 764704289: 211 dropped running frames (~0.1%) (in pre-processing).

Finished creating session 764704289.

Creating session 765193831...

Loading stimulus and alignment info...

Loading ROI trace info...

WARNING: Session 765193831: 4 noisy ROIs (mean below 0 or median above midrange) are also included in the NaN ROI attributes (but not set to NaN): 3, 63, 88, 134.

Loading running info...

WARNING: Session 765193831: 345 dropped running frames (~0.1%) (in pre-processing).

Finished creating session 765193831.

Creating session 766502238...

Loading stimulus and alignment info...

Loading ROI trace info...

WARNING: Session 766502238: 4 noisy ROIs (mean below 0 or median above midrange) are also included in the NaN ROI attributes (but not set to NaN): 18, 114, 136, 240.

Loading running info...

WARNING: Session 766502238: 387 dropped running frames (~0.2%) (in pre-processing).

Finished creating session 766502238.

Creating session 777914830...

Loading stimulus and alignment info...

Loading ROI trace info...

WARNING: Session 777914830: 1 noisy ROIs (mean below 0 or median above midrange) are also included in the NaN ROI attributes (but not set to NaN): 45.

Loading running info...

WARNING: Session 777914830: 381 dropped running frames (~0.2%) (in pre-processing).

Finished creating session 777914830.

Creating session 778864809...

Loading stimulus and alignment info...

Loading ROI trace info...

Loading running info...

WARNING: Session 778864809: 630 dropped running frames (~0.3%) (in pre-processing).

Finished creating session 778864809.

Creating session 758519303...

Loading stimulus and alignment info...

Loading ROI trace info...



Loading running info...

WARNING: Session 758519303: 175 dropped running frames (~0.1%) (in pre-processing).

Finished creating session 758519303.

Then run through the sessions and do whatever with them.

Session ID: 764704289 (mouse 6, session 1)

visflow: 33 sequences

gabors: 96 sequences

Session ID: 765193831 (mouse 6, session 2)

visflow: 34 sequences

gabors: 98 sequences

Session ID: 766502238 (mouse 6, session 3)

visflow: 29 sequences

gabors: 94 sequences

Session ID: 777914830 (mouse 8, session 1)

visflow: 32 sequences

gabors: 83 sequences

Session ID: 778864809 (mouse 8, session 2)

visflow: 29 sequences

gabors: 88 sequences

Session ID: 758519303 (mouse 1, session 1)

visflow: 31 sequences

gabors: 94 sequences

#### 1.4.2 5. Retrieving ROI masks from session.

Boolean ROI masks can be obtained for Session.

For **dendritic sessions**, the Session is built to assume that **EXTRACT** (not **allen**) ROI data is to be used. This can be checked by checking `self.dend`. As long as `self.dend` is properly set, the correct masks will be loaded.

Dendritic session, ROI type: `extr`


Somatic session, ROI type: `allen`

Masks can be loaded as follows, with dimensions: **ROI x height x width**, retrieving only masks for ROIs that are valid for dF/F traces.


Or if all masks are needed, the attributes can simply be retrieved.

One way to check which ROIs are not valid, is using `self.get_nanrois()`

This is a tool to visualize ROIs, where specific ROIs can be set to red using a `valid_mask`.



This is a tool to visualize ROI contours, optionally localized around an ROI of interest.



### 1.4.3 6. Last notes

List the methods/attributes attached to `Session` and `Stim` objects.

Session: `Session (758519303)`

`_dend`, `_extract_sess_attribs`, `_full_table`, `_get_roi_facts`, `_init_directory`,  
`_init_roi_facts_df`, `_load_stim_df`, `_load_stim_dict`, `_load_stims`, `_nanrois_dff`,  
`_nrois`, `_nwb`, `_only_tracked_rois`, `_roi_masks`, `_set_dend_type`, `_set_nanrois`,  
`_set_nanrois_tracked`, `_set_roi_attributes`, `_set_tracked_rois`, `_stim2twopfr`,

align\_pkl, all\_files, any\_files, behav\_video\_h5, check\_flanks, convert\_frames, correct\_data\_h5, data\_loaded, date, dend, depth, dir, drop\_tol, expdir, expid, extract\_info, gabors, get\_active\_rois, get\_fr\_ran, get\_frames\_timestamps, get\_nanrois, get\_nrois, get\_plateau\_roi\_traces, get\_pup\_data, get\_roi\_masks, get\_roi\_seqs, get\_roi\_traces, get\_run\_velocity, get\_run\_velocity\_by\_fr, get\_single\_roi\_trace, get\_stim, grayscr, home, line, load\_pup\_data, load\_roi\_info, load\_run\_data, mouse\_df, mouse\_dir, mouse\_n, mouseid, n\_stims, notes, nwb, only\_tracked\_rois, pass\_fail, plane, procdir, pup\_data\_h5, pup\_video\_h5, roi\_extract\_json, roi\_facts\_df, roi\_mask\_file, roi\_masks, roi\_names, roi\_objectlist, roi\_trace\_dff\_h5, roi\_trace\_h5, run\_data, runtime, segid, sess\_n, sessid, set\_only\_tracked\_rois, stim2twopfr, stim\_df, stim\_fps, stim\_pkl, stim\_seed, stim\_sync\_h5, stims, stimtypes, time\_sync\_h5, tot\_stim\_fr, tot\_twop\_fr, tracked\_rois, twop2stimfr, twop\_fps, visflow, zstack\_h5

Gabors: Session (758519303)

\_set\_block\_params, all\_gabfr, all\_gabfr\_mean\_oris, block\_params, deg\_per\_pix, exp\_gabfr, exp\_gabfr\_mean\_oris, exp\_max\_s, exp\_min\_s, get\_A\_frame\_1s, get\_A\_segs, get\_all\_unexp\_segs, get\_all\_unexp\_stim\_fr, get\_fr\_by\_seg, get\_frames\_by\_criteria, get\_n\_fr\_by\_seg, get\_pup\_diam\_data, get\_pup\_diam\_stats\_df, get\_roi\_data, get\_roi\_stats\_df, get\_run, get\_run\_data, get\_run\_stats\_df, get\_segs\_by\_criteria, get\_segs\_by\_frame, get\_start\_unexp\_segs, get\_start\_unexp\_stim\_fr\_trans, get\_stats\_df, get\_stim\_beh\_sub\_df, get\_stim\_df\_by\_criteria, get\_stim\_par\_by\_frame, get\_stim\_par\_by\_seg, kappas, n\_patches, n\_segs\_per\_seq, ori\_ran, phase, seg\_len\_s, sess, sf, size\_ran, stim\_fps, stimtype, unexp\_gabfr, unexp\_gabfr\_mean\_oris, unexp\_max\_s, unexp\_min\_s, win\_size

Visflow: Session (758519303)

\_set\_block\_params, block\_params, deg\_per\_pix, exp\_max\_s, exp\_min\_s, get\_all\_unexp\_segs, get\_all\_unexp\_stim\_fr, get\_dir\_segs\_exp, get\_fr\_by\_seg, get\_frames\_by\_criteria, get\_n\_fr\_by\_seg, get\_pup\_diam\_data, get\_pup\_diam\_stats\_df, get\_roi\_data, get\_roi\_stats\_df, get\_run, get\_run\_data, get\_run\_stats\_df, get\_segs\_by\_criteria, get\_segs\_by\_frame, get\_start\_unexp\_segs, get\_start\_unexp\_stim\_fr\_trans, get\_stats\_df, get\_stim\_beh\_sub\_df, get\_stim\_df\_by\_criteria, main\_flow\_dirs, n\_squares, prop\_flipped, seg\_len\_s, sess, speed, square\_sizes, stim\_fps, stimtype, unexp\_max\_s, unexp\_min\_s, win\_size

Grayscr: Session (758519303)

get\_all\_fr, get\_start\_fr, get\_stop\_fr, sess

List Stim object attribute values.

Gabors: Gabors (stimulus from session 758519303)

```
{'all_gabfr': ['A', 'B', 'C', 'D', 'U'],
 'all_gabfr_mean_oris': [0.0, 45.0, 90.0, 135.0, 180.0, 225.0],
 'block_params':
   start_seg  stop_seg  num_segs  start_time_sec  stop_time_sec
duration_sec \
0           1       6801       6800       43.540279    2085.237221    2041.696942
```

	start_frame_stim	stop_frame_stim	num_frames_stim	start_frame_twop \
0	1800	124200	122400	1045

	stop_frame_twop	num_frames_twop	gabor_kappa
0	62456	61411	16 ,

```

'deg_per_pix': 0.06251912565744862,
'exp_gabfr': ['A', 'B', 'C', 'D'],
'exp_gabfr_mean_oris': [0.0, 45.0, 90.0, 135.0],
'exp_max_s': 90,
'exp_min_s': 30,
'kappas': [16],
'n_patches': 30,
'n_segs_per_seq': 5,
'ori_ran': [0, 360],
'phase': 0.25,
'seg_len_s': 0.3,
'sess': Session (758519303),
'sf': 0.04,
'size_ran': [159.9510532951381, 319.9021065902762],
'stim_fps': 59.95049782968851,
'stimtype': 'gabors',
'unexp_gabfr': ['U'],
'unexp_gabfr_mean_oris': [90.0, 135.0, 180.0, 225.0],
'unexp_max_s': 6,
'unexp_min_s': 3,
'win_size': [1920, 1200]}

```

Visflow: Visflow (stimulus from session 758519303)

	start_seg	stop_seg	num_segs	start_time_sec	stop_time_sec
0	6802	7822	1020	2115.262241	3136.110155
1	7823	8843	1020	3166.135568	4186.983399

	start_frame_stim	stop_frame_stim	num_frames_stim	start_frame_twop \
0	126000	187200	61200	63359
1	189000	250200	61200	94967

	stop_frame_twop	num_frames_twop	main_flow_direction	square_size \
0	94064	30705	right (temp)	128
1	125671	30704	left (nasal)	128

	square_number
0	105
1	105 ,

```

'deg_per_pix': 0.06251912565744862,
'exp_max_s': 90,
'exp_min_s': 30,
'main_flow_direcs': ['left (nasal)', 'right (temp)'],

```

```
'n_squares': [105],
'prop_flipped': 0.25,
'seg_len_s': 1,
'sess': Session (758519303),
'speed': 799.7552664756905,
'square_sizes': [128],
'stim_fps': 59.95049782968851,
'stimtype': 'visflow',
'unexp_max_s': 4,
'unexp_min_s': 2,
'win_size': [1920, 1200]}

Grayscr: Grayscr (session 758519303)
{'sess': Session (758519303)}
```